



Test kit, phosphate, model PO-23

Product #:

224902

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	F
Measuring range:	Dual Range
Method:	Color disc/Ascorbic acid for clear samples
Model:	PO-23
Number of tests:	100
Parameter:	Phosphate
Platform:	Color Disc
Range:	0.1 - 4 mg/L PO ₄
Range 2:	1 - 40 mg/L PO ₄
Smallest increment:	0.1
	1
Standard range:	0.1 - 4 mg/L PO ₄ ; 1 - 40 mg/L PO ₄
Test Kit Model:	PO-23