



Test kit, phosphate, total ortho-/meta-, model PO-24

Product #:

225001

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Kit provides measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

Kit provides rapid measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

Specifications

Case Style:	F
Measuring range:	Triple Range
Method:	Color disc/Ascorbic acid
Model:	PO-24
Number of tests:	100 Clear
	100 Turbid
Parameter:	Phosphate
Platform:	Color Disc
Range:	0 - 0.8 mg/L PO_4
Range 2:	0 - 4 mg/L PO_4
Range 3:	0 - 40 mg/L PO_4
Smallest increment:	0.02
	0.1
	1
Standard range:	0 - 0.8 mg/L PO_4 ; 0 - 4 mg/L PO_4 ; 0 - 40 mg/L PO_4
Test Kit Model:	PO-24