



Test Kit Iron, 0 - 4 mg/L Fe, Colour disc method, 100 tests

Product #: 146400
ZAR Price: Contact Hach
Ships within 3 days

Specifications

Case Style:	B
Measuring range:	Medium Single Range
Method:	Color disc/ 1,10 Phenanthroline
Model:	IR-18
Number of tests:	100
Parameter:	Iron
Platform:	Color Disc
Range:	0 - 4 mg/L Fe
Smallest increment:	0.1
Standard range:	0 - 4 mg/L Fe
Test Kit Model:	IR-18
Type:	Single Metric