



Iron color disc test kit, model IR-18B

Product #: 146401
ZAR Price: Contact Hach
Ships within 3 days

Using the simple color disc method, this kit measures Iron (medium range - 0.2-7 mg/L) and contains reagent for 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	B
Measuring range:	High Single Range
Method:	Color disc/ 1,10 Phenanthroline
Model:	IR-18B
Number of tests:	100
Parameter:	Iron
Platform:	Color Disc
Range:	0.25 - 7 mg/L Fe
Smallest increment:	0.25
Standard range:	0.25 - 7 mg/L Fe
Test Kit Model:	IR-18B
Type:	Single Metric