



## Iron color disc test kit, model IR-18A

**Product #:** 146500  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

Using the simple color disc method, this kit measures Iron (Range - 0.02-1 mg/L) and contains 100 tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Case Style:	B
Measuring range:	Low Single Range
Method:	Color disc/ 1,10 Phenanthroline
Model:	IR-18A
Number of tests:	100
Parameter:	Iron
Platform:	Color Disc
Range:	0.02 - 1.00 mg/L Fe
Smallest increment:	0.02
Standard range:	0.02 - 1.00 mg/L Fe
Test Kit Model:	IR-18A
Type:	Single Metric