



## Chromium color disc test kit, model CH-8

**Product #:** 183400  
**ZAR Price:** Contact Hach  
**Available**

Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the color disc method to determine chromium (Range: 0.1-1.5 mg/L), this kit contains 100 tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Case Style:	D
Method:	Color disc/  Diphenylcarbazide
Model:	CH-8
Number of tests:	100
Parameter:	Chromium Hexavalent
Platform:	Color Disc
Range:	0 - 1.3 Cr <sup>6+</sup>
Smallest increment:	0.1
Test Kit Model:	CH-8