



Detergents test kit, model DE-2

Product #:

143203

ZAR Price:

Contact Hach

Ships within 1 week

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The presence of detergents in surface or ground waters indicates industrial or domestic pollution. In natural waters anionic LAS and ABS detergents cause excessive foaming, produce undesirable taste and odors, promote growth of disease-producing organisms, and tend to kill aquatic life. This kit's simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	Q
Method:	Color disc/ Toluidine blue-O Chloroform
Model:	DE-2
Number of tests:	32
Parameter:	Detergents - As LAS and/or ← ABS detergents
Platform:	Color Disc
Range:	0 - 1.2 mg/L LAS
Test Kit Model:	DE-2