



## Hydrazine color disc test kit, model HY-2

**Product #:**

**184900**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method with a range of 0.02-1 mg/L, this kit contains 300 tests.

**Measures Hydrazine, a common oxygen scavenger in boiler water, using the simple color disc method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Case Style:	D
Method:	Color disc/  p-Dimethylamino-  benzaldehyde
Model:	HY-2
Number of tests:	300
Parameter:	Detergents
Platform:	Color Disc
Range:	0 - 1 mg/L $N_2H_4$
Smallest increment:	0
Standard range:	0 - 1 mg/L $N_2H_4$
Test Kit Model:	HY-2