



## Test kit, manganese color disc, model MN-PAN

**Product #:**

**2350800**

**ZAR Price:**

**Contact Hach**

**Ships within 3-5 weeks**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Case Style:	D
Method:	Color disc/PAN
Model:	MN-PAN
Number of tests:	50
Parameter:	Manganese
Platform:	Color Disc
Range:	0 - 0.70 mg/L Mn
Smallest increment:	0.05
Standard range:	0 - 0.70 mg/L Mn
Test Kit Model:	MN-PAN