



Test kit cyanide, model CYN-3

Product #:	201002
ZAR Price:	Contact Hach
Ships within 2 weeks	

The main sources of cyanide, an extremely poisonous substance, are metal-cleaning and electroplating baths, gas scrubber, gas works, coke ovens, and various chemical processes. It can be removed easily from waste effluents through alkaline or neutral chlorination procedures.

Simple, low-cost method using a stepped color comparator

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	R
Method:	Color disc/ Pyridine-Pyrazolone
Model:	CYN-3
Number of tests:	100
Parameter:	Cyanide
Platform:	Color Disc
Range:	0 - 0.3 mg/L CN ⁻
Smallest increment:	0.02
Standard range:	0 - 0.3 mg/L CN ⁻
Test Kit Model:	CYN-3