



Test kit, free chlorine, model CN-70F

Product #: 1454201
ZAR Price: Contact Hach
Call for ship date

Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Quickly and accurately measure free chlorine levels in the field or in the lab

Specifications

Case Style:	D
Measuring range:	Dual Range
Method:	Color disc/DPD
Number of tests:	200
Parameter:	Chlorine Free
Platform:	Color Disc
Range:	0.02 - 0.68 mg/L Cl ₂
Range 2:	0 - 3.4 mg/L Cl ₂
Smallest increment:	0.02
	0.1
Test Kit Model:	CN-70F ²