



Test Kit, chlorine (free and total), model CN-80

Product #: 2129000
ZAR Price: Contact Hach
Ships within 3 days

Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	D
Measuring range:	Triple Range
Method:	Color Disc/DPD
Model:	CN-80
Number of tests:	100
Parameter:	Chlorine (Free and Total)
Platform:	Color Disc
Range:	0 - 0.68 mg/L Cl_2
Range 2:	0 - 3.4 mg/L Cl_2
Range 3:	0 - 10 mg/L Cl_2
Smallest increment:	0.02
	0.1
	0.5
Standard range:	0 - 0.68 mg/L Cl_2 ; 0 - 3.4 mg/L Cl_2 ; 0 - 10 mg/L Cl_2
Test Kit Model:	CN-80 ²