



Chlorine (free & total) colour disc test kit, model CN-66 (pk/100)

Product #:

223101

ZAR Price:

Contact Hach

Ships within 3 days

Measuring range 0.1-3.4 mg/L Cl₂

Specifications

Case Style:	B
Measuring range:	Single Range
Method:	Color disc/DPD
Model:	CN-66
Number of tests:	50 Free, 50 Total
Parameter:	Chlorine (Free and Total)
Platform:	Color Disc
Range:	0 - 3.4 mg/L Cl ₂
Smallest increment:	0.2
Standard range:	0 - 3.4 mg/L Cl ₂
Test Kit Model:	CN-66 ²