



Chlorine, free; colour disc test kit, model CN-66F (pk/100)

Product #: 223102
ZAR Price: Contact Hach
Ships within 3 days

Measuring range 0.1 - 3.4 mg/L Cl₂

Specifications

Case Style:	B
Footnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Measuring range:	Single Range
Method:	Color disc/DPD
Model:	CN-66F
Number of tests:	100
Parameter:	Chlorine Free
Platform:	Color Disc
Range:	0 - 3.4 mg/L Cl ₂
Smallest increment:	0.2
Standard range:	0 - 3.4 mg/L Cl ₂
Test Kit Model:	CN-66F ²