



Silica Test Kit, SI-7, 0.05 - 1.0 mg/L SiO₂

Product #:

2255000

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Simple, low-cost colour disc method for Silica determination.

Specifications

Case Style:	D
Measuring range:	Single Range
Model:	SI-7
Number of tests:	100
Parameter:	Silica
Platform:	Color Disc
Range:	0 - 1 mg/L SiO ₂
Smallest increment:	0.05 g/L
Standard range:	0 - 1 mg/L SiO ₂
Test Kit Model:	SI-7