



AccuVac kit, dissolved oxygen

Product #: 2515050
ZAR Price: Contact Hach
Ships within 3 days

Low DO levels usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. The simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	G
Method:	Color disc/HRDO
Model:	Color disc/HRDO
Number of tests:	25
Parameter:	Dissolved Oxygen - As DO
Platform:	Accuvac
Range:	0 - 15 mg/L O ₂
Smallest increment:	0.2
Test Kit Model:	AccuVac® Kit