



Test kit, ozone, model OZ-2

Product #: 2064400
ZAR Price: Contact Hach
Ships within 3 days

Hach's most economical test kit for ozone measurement

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	B
Method:	Color disc/DPD
Model:	OZ-2
Number of tests:	100
Parameter:	Ozone
Platform:	Color Disc
Range:	0 - 2.2 mg/L O ₃
Smallest increment:	0.05
Standard range:	0 - 2.2 mg/L O ₃
Test Kit Model:	OZ-2