



## Phosphonate Test Kit, Model PN-10, with 115 Vac UV Lamp and Power Supply

**Product #:**

**2113300**

**ZAR Price:**

Contact Hach

**Ships within 1 week**



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Contains 100 tests.

**Simple, fast, color disc method**

**UV lamp digestion converts phosphonates into reactive phosphate**

**Control phosphonate concentration to save chemical costs**

**Hach's reliable PhosVer 3 reagent provides results in under 15 minutes**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Case Style:	F
Method:	Color disc/UV photolysis (115 Vac, 60 Hz)
Model:	PN-10
Number of tests:	0 - 100
Parameter:	Phosphonate - As $\text{PO}_4^{3-}$
Platform:	Color Disc
Power supply:	115 Vac UV Lamp and Power Supply
Range:	0 - 100 mg/L $\text{PO}_4$
Smallest increment:	1
Test Kit Model:	PN-10
What's included?:	Includes color disc, comparator box and tubes, volumetric flask, graduated cylinder, reaction bottle, reagents, UV lamp, 115 V power supply, UV safety goggles, instructions and carrying case.

---

### What's included?

Includes color disc, comparator box and tubes, volumetric flask, graduated cylinder, reaction bottle, reagents, UV lamp, 115 V power supply, UV safety goggles, instructions and carrying case.