



Test kit, phosphonate, model PN-10, w/out uv lamp & power supply

Product #:

2113301

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Contains 100 tests.

Simple, fast, color disc method

UV lamp digestion converts phosphonates into reactive phosphate

Control phosphonate concentration to save chemical costs

Hach's reliable PhosVer 3 reagent provides results in under 15 minutes

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	F
Method:	Color disc
Model:	PN-10
Number of tests:	0 - 100
Parameter:	Phosphonate
Platform:	Color Disc
Power supply:	No UV Lamp Nor Power Supply
Range:	0 - 100 mg/L PO_4
Smallest increment:	1
Standard range:	0 - 100 mg/L PO_4
Test Kit Model:	PN-10

(As above but without UV lamp and power supply.)