



Triazole Test Kit, Model TZ-1, with 115 Vac UV Lamp and Power Supply

Product #: 2167500
ZAR Price: Contact Hach
Ships within 1 week

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

Measures common corrosion and biological growth inhibitor using simple color disc method

UV lamp digestion

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	F
Method:	Color disc/ UV photolysis or tolyltriazole (115 Vac, 60 Hz)
Model:	TZ-1
Number of tests:	50
Parameter:	Triazole - As benzotriazole
Platform:	Color Disc
Power supply:	115 Vac UV Lamp and Power Supply
Range:	0 - 15.0 mg/L
Smallest increment:	0.5
Test Kit Model:	TZ-1
What's included?:	Includes benzotriazole and tolyltriazole color discs, comparator box and tubes, reaction bottle, reagents, UV lamp, 115 V power supply, UV safety goggles, instructions and carrying case.

What's included?

Includes benzotriazole and tolyltriazole color discs, comparator box and tubes, reaction bottle, reagents, UV lamp, 115 V power supply, UV safety goggles, instructions and carrying case.