



Test kit, triazole, model TZ-1, w/ 230 vac uv lamp & power supply

Product #: 2167502
ZAR Price: Contact Hach
Ships within 2 weeks

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

Measures common corrosion and biological growth inhibitor using simple color disc method

UV lamp digestion

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	F
Method:	Color disc/ UV photolysis (230 Vac, 50 Hz)
Model:	TZ-1
Number of tests:	50
Parameter:	Triazole -As benzotriazole
Platform:	Color Disc
Power supply:	230 Vac UV Lamp and Power Supply
Range:	0 - 15.0 mg/L
Smallest increment:	0.5
Test Kit Model:	TZ-1