



Iron Test Kit, Low Range (Model IR-24)

Product #:

255600

ZAR Price:

Contact Hach

Ships within 2 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Simple, low-cost colour disc method for Iron determination. The Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

Specifications

Case Style:	D
Measuring range:	Dual Range
Method:	Color disc/FerroZine
Model:	IR-24
Number of tests:	100
Parameter:	Iron Ferrous
Platform:	Color Disc
Range:	0 - 0.2 mg/L Fe
Range 2:	0 - 1.0 mg/L Fe
Smallest increment:	0.002
	0.01
Standard range:	0 - 0.2 mg/L Fe; 0 - 1.0 mg/L Fe
Test Kit Model:	IR-24
Type:	Double Metric