



## Iron (Ferrous) Color Disc Test Kit, Model IR-18C, 0-7 mg/L

**Product #:** 2667200  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

The most accurate visual test kit method, Hach color discs provide fast, accurate comparisons. Range: 0-7 mg/L. Contains 100 tests.

**Simple, fast, color disc method**

**Rugged test kit contains all required reagents and apparatus to perform analysis in the field**

---

### Specifications

Case Style:	B
Measuring range:	High Single Range
Method:	Color disc/ 1,10 Phenanthroline
Model:	IR-18C
Number of tests:	100
Parameter:	Iron Ferrous
Platform:	Color Disc
Range:	0 - 7 mg/L Fe <sup>2+</sup>
Smallest increment:	0.2
Test Kit Model:	IR-18C
Type:	Single Metric