



## Manganese test kit, model MN-5

**Product #:**

**146700**

ZAR Price:  
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. An economical solution, this kit uses the simple, accurate color disc method.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Specifications

Case Style:	D
Method:	Color disc/Cold periodate
Model:	MN-5
Number of tests:	100
Parameter:	Manganese
Platform:	Color Disc
Range:	0 - 3.0 mg/L Mn
Smallest increment:	0.1
Standard range:	0 - 3.0 mg/L Mn
Test Kit Model:	MN-5