



Test kit, oxygen scavenger, color disc, model OS-1, 100 tests

Product #:

2349500

ZAR Price:

Contact Hach

Ships within 3-5 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 100 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	R
Method:	Color disc/Iron reduction
Model:	OS-1
Number of tests:	50
Parameter:	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Platform:	Color Disc
Range:	0 - 52 µg/L DEHA (and other Oxygen Scavengers)
Range 2:	0 - 1360 µg/L DEHA (and other Oxygen Scavengers)
Range 3:	0 - 300 µg/L DEHA (and other Oxygen Scavengers)
Test Kit Model:	OS-1