



Test kit, phenols, model PL-1

Product #: 2483600
ZAR Price: Contact Hach
Ships within 3-5 weeks



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Utilizing the simple color disc method with a range of 0 - 4 mg/L, this kit contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	D
Method:	Color disc/ 4-aminoantipyrine
Model:	PL-1
Number of tests:	100
Parameter:	Phenols
Platform:	Color Disc
Range:	0 - 4.0 mg/L
Smallest increment:	0.2 mg/L
Standard range:	0 - 4.0 mg/L
Test Kit Model:	PL-1