



Test kit, glutaraldehyde, model GT-1

Product #:

2587200

ZAR Price:

Contact Hach

Ships within 3-5 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

From oil fields to cooling towers, hospitals to food plants, poultry farms to dental offices, the Glutaraldehyde Test Kit can help monitor disinfection processes. The test kit measures from 0.5 to 4000 ppm glutaraldehyde.

Helps monitor disinfection processes

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	Q
Method:	Color disc
Model:	GT-1
Number of tests:	100
Parameter:	Glutaraldehyde
Platform:	Color Disc
Range:	1 - 2800 mg/L multiple ranges
Smallest increment:	0.5
Standard range:	1 - 2800 mg/L multiple ranges
Test Kit Model:	GT-1