



## Phosphorus, orthophosphate (reactive) test kit, model PO-14

**Product #:**

**147500**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Simple, low-cost color disc method is the most accurate visual test kit method

---

### Specifications

Case Style:	D
Measuring range:	Single Range
Method:	Color disc/Stannous
Model:	PO-14
Number of tests:	100 C**
Parameter:	Phosphorus
Platform:	Color Disc
Range:	0 - 4.4 mg/L PO <sub>4</sub>
Smallest increment:	0.1
	1
Standard range:	0 - 4.4 mg/L PO <sub>4</sub>
Test Kit Model:	PO-14