



## Total chlorine test kit (colour disk)

**Product #:** 1454202  
**ZAR Price:** Contact Hach  
**Available**

Colour comparison test for easy and low-cost determination of total chlorine. Measuring range 0.02 - 0.7 / 0 - 3.4 mg/L Cl<sub>2</sub>.

---

### Specifications

|                     |                                                 |
|---------------------|-------------------------------------------------|
| Case Style:         | D                                               |
| Measuring range:    | Dual Range                                      |
| Method:             | Color disc/DPD                                  |
| Number of tests:    | 100 Total                                       |
| Parameter:          | Chlorine Total                                  |
| Platform:           | Color Disc                                      |
| Range:              | 0.02 - 3.4 mg/L Cl <sub>2</sub> Cl <sub>2</sub> |
| Range 2:            | 0.02 - 0.68 mg/L Cl <sub>2</sub>                |
| Smallest increment: | 0.02                                            |
|                     | 0.1                                             |
| Test Kit Model:     | CN-70T <sup>2</sup>                             |