



## Test kit color disc chlorine dioxide, model CLD-2

**Product #:** 2842800  
**ZAR Price:** Contact Hach  
**Ships within 2 weeks**

Monitor  $\text{ClO}_2$  levels on-site with the NEW Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable and EPA-accepted DPD method to monitor  $\text{ClO}_2$  concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

**Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant**

**USEPA-accepted for waste water and potable water**

**Approved by some states for potable water under specific conditions (check with your regulatory agency)**

---

### Specifications

Case Style:	D
Method:	Color Disc/DPD
Model:	CLD-2
Number of tests:	100
Parameter:	Chlorine Dioxide
Platform:	Color Disc
Range:	0 - 2.0 ppm $\text{ClO}_2$
Smallest increment:	0.1
Standard range:	0 - 2.0 ppm $\text{ClO}_2$
Test Kit Model:	CLD-2 <sup>2</sup>