



Test kit, nitrate-nitrite, model NI-12

Product #: 1408100
ZAR Price: Contact Hach
Ships within 3-5 weeks

Simple, low-cost color disc method is the most accurate visual test kit method.

Specifications

Case Style:	D
Measuring range:	Dual Range
Method:	Color disc/Cadmium reduction, Diazotization
Model:	NI-12
Number of tests:	100
Parameter:	Nitrogen, Nitrate/Nitrite
Platform:	Color Disc
Range:	0 - 40 mg/L NO_3^- -N
Range 2:	0.01 - 0.40 mg/L NO_2^- -N
Smallest increment:	1 0.01
Standard range:	0 - 40 mg/L NO_3^- -N; 0.01 - 0.40 mg/L NO_2^- -N
Test Kit Model:	NI-12**