



Test kit, nitrogen, nitrate color disc, model NI-14

Product #: 1416100
ZAR Price: Contact Hach
Ships within 3 days

Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO₃ and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	B
Measuring range:	Dual Range
Method:	Color disc/ Cadmium reduction
Model:	NI-14
Number of tests:	50
Parameter:	Nitrogen, Nitrate
Platform:	Color Disc
Range:	0 - 1 mg/L NO ₃ ⁻ -N
Range 2:	0 - 10 mg/L NO ₃ ⁻ -N
Smallest increment:	0.02 0.2
Standard range:	0 - 1 mg/L NO ₃ ⁻ -N; 0 - 10 mg/L NO ₃ ⁻ -N
Test Kit Model:	NI-14**