



## Nitrate test kit, model NI-11

**Product #:** 146803  
**ZAR Price:** Contact Hach  
**Available**

One of Hach's most popular test kits, it uses the simple, low-cost color disc method.

**One of Hach's most popular test kits**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Specifications

Case Style:	B
Measuring range:	Single Range
Method:	Color disc/ Cadmium reduction
Model:	NI-11
Number of tests:	100
Parameter:	Nitrogen, Nitrate
Platform:	Color Disc
Range:	0 - 40 mg/L $\text{NO}_3\text{-N}$
Smallest increment:	1.0
Standard range:	0 - 40 mg/L $\text{NO}_3\text{-N}$
Test Kit Model:	NI-11**