



Test kit, nitrite, nitrogen, model NI-15

Product #:	2182000
ZAR Price:	Contact Hach
Ships within 3 days	

Simple, low-cost color disc method is the most accurate visual test kit method.

Specifications

Case Style:	D
Measuring range:	Single Range
Method:	Color disc/Diazotization
Model:	NI-11
Number of tests:	100
Parameter:	Nitrite
Platform:	Color Disc
Range:	0 - 0.40 mg/L NO_2^- -N
Smallest increment:	0.01
Standard range:	0 - 0.40 mg/L NO_2^- -N
Test Kit Model:	NI-15