



## Test kit, nitrogen , model NI-8

**Product #:** 224100  
**ZAR Price:** Contact Hach  
**Available**



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

**Contains mercury. Dispose of according to local, state and federal regulations.**

Designed for testing ammonia in freshwater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

**One of Hach's most popular test kits**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Specifications

Case Style:	B
Method:	Color disc/Nessler reagent
Number of tests:	100
Parameter:	Nitrogen, Ammonia, mid range - For freshwater
Platform:	Color Disc
Range:	0 - 3.0 mg/L NH <sub>3</sub> -N
Smallest increment:	0.1
Test Kit Model:	NI-8