



Test kit, nitrogen, ammonia, model NI-SA

Product #:

2428700

ZAR Price:

Contact Hach

Ships within 1 week

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Designed for testing ammonia in freshwater or seawater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

Designed for testing ammonia in freshwater or seawater

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	D
Method:	Color disc/Salicylate
Number of tests:	100
Parameter:	Nitrogen, Ammonia, mid range - For freshwater or seawater**
Platform:	Color Disc
Range:	0 - 2.4 mg/L NH ₃ -N
Smallest increment:	0.1
Test Kit Model:	NI - SA