



pH test kit, 5.5 - 8.5 pH, model 17F

Product #:

147006

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	B
Method:	Color disc/ Bromthymol blue
Model:	17F
Number of tests:	200
Parameter:	pH, mid range
Platform:	Color Disc
Range:	5.6 - 8.4 pH
Smallest increment:	0.1 pH
Test Kit Model:	17F