



pH test kit, 6.5 - 8.5 pH, model 17H

Product #:	147008
ZAR Price:	Contact Hach
Ships within 3 days	

Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	B
Method:	Color disc/Phenol red
Model:	17H
Number of tests:	200
Parameter:	pH, mid range
Platform:	Color Disc
Range:	6.6 - 8.4 pH
Smallest increment:	0.1 pH
Test Kit Model:	17H