



Test kit, phosphonate, model PN-10, w/ 230 vac uv lamp & power supply

Product #:

2113302

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Contains 100 tests.

Simple, fast, color disc method

UV lamp digestion converts phosphonates into reactive phosphate

Control phosphonate concentration to save chemical costs

Hach's reliable PhosVer 3 reagent provides results in under 15 minutes

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

| | |
|---------------------|--|
| Case Style: | F |
| Method: | Color disc/UV photolysis (230 Vac, 50 Hz) |
| Model: | PN-10 |
| Number of tests: | 0 - 100 |
| Parameter: | Phosphonate - As PO_4^{3-} |
| Platform: | Color Disc |
| Power supply: | 230 Vac UV Lamp and Power Supply |
| Range: | 0 - 100 mg/L PO_4 |
| Smallest increment: | 1 |
| Test Kit Model: | PN-10 |