



Test Kit, Accuvac Color Disc, Ozone, 0.05 - 1.50 mg/L

Product #: 2518050
ZAR Price: Contact Hach
Ships within 2 weeks

This ozone test kit uses the fast, accurate visual method for determining ozone concentrations. Use to control disinfectant levels in water, wastewater, and industrial water. Kit contains 25 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

AccuVac vials mean no mess, no measuring

Indigo trisulfonate method masks interferences from chlorine and other oxidizing agents and makes ozone analysis easier

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Specifications

Case Style:	G
Method:	Color disc/ Indigo trisulfonate
Model:	NA
Number of tests:	25
Parameter:	Ozone
Platform:	Accuvac
Range:	0 - 1.50 mg/L O ₃
Smallest increment:	0.05
Standard range:	0 - 1.50 mg/L O ₃
Test Kit Model:	AccuVac® Kit