



Test kit, copper (free & total) color disc, model CU-6

Product #:

2194100

ZAR Price:

Contact Hach

Ships within 2 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Copper may exist in natural waters and effluents as a soluble salt or as suspended solids. A small amount is essential for plants and animals, and concentrations exceeding 0.1 mg/L are useful for controlling algae and plankton growth. Quantities ranging from 0.02 to 0.1 mg/L are toxic for some fish, so its use for treating fish disease requires careful monitoring. Range: 0.1-4 mg/L; Contains 50 Free and 50 Total tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	B
Method:	Color-disc/Bicinchoninate, Hydrosulfite reduction
Model:	CU-6
Number of tests:	50 (Free), 50 (Total)
Parameter:	Copper (free & total)
Platform:	Color Disc
Range:	0 - 4 mg/L Cu
Smallest increment:	0.1
Test Kit Model:	CU-6