



Iron Color Disc Test Kit, Model IR-10

Product #: 146430
ZAR Price: Contact Hach
Ships within 3-5 weeks

Using the simple color disc method, this kit measures Iron (medium range: 0.1 - 10 mg/L) and contains reagent for 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Quickly and accurately measure iron in the field or in the lab

Specifications

Case Style:	B
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Method:	Color disc/ 1,10 Phenanthroline
Number of tests:	100
Parameter:	Iron, medium range - as Fe
Platform:	Color Disc
Range:	0 - 10 mg/L
Test Kit Model:	IR-10
What's included?:	Includes one reagent, two color discs, comparator box, two viewing tubes, instructions and carrying case.

What's included?

Includes one reagent, two color discs, comparator box, two viewing tubes, instructions and carrying case.