



Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

Product #:

143601

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO_2 , with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Specifications

Case Style:	D
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Measuring range:	2 - 40 mg/L
Method:	Drop count titration/ Sodium hydroxide Phenolphthalein
Number of tests:	200
Parameter:	Carbon Dioxide - As CO_2
Platform:	Drop Count Titration
Range 2:	1.25 - 25 mg/L
Range 3:	5 - 100 mg/L CO_2
Smallest increment:	1.25 2 5
Test Kit Model:	CA-23