



Chloride Low Range Test Kit, Model 8-P

Product #:

144001

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Using silver nitrate titrant, this kit provides fast, accurate detection of chloride.

Low cost and easy to use drop count titration method using a silver nitrate titrant

Ideal for quickly determining quantitative results

Specifications

Case Style:	B
Measuring range:	5 - 100 mg/L & 20 - 400 mg/L Cl^-
Method:	Drop count titration/ Silver nitrate
Number of tests:	100
Parameter:	Chloride - As Cl^-
Platform:	Drop Count Titration
Smallest increment:	5
Test Kit Model:	8-P