



TEST KIT, 5-EP mg/L HARDNESS

Product #: 145401
ZAR Price: Contact Hach
Ships within 3 days

Total Hardness drop count titration/EDTA test kit with range of 20 - 400 mg/L; contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	B
Measuring range:	20 - 400 mg/L CaCO_3
Method:	Drop count titration/ EDTA-powder pillows
Number of tests:	100
Parameter:	Hardness, total - As CaCO_3
Platform:	Drop Count Titration
Smallest increment:	20
Test Kit Model:	5-EP MG-L