



Hardness (Total & Calcium) Test Kit, Model HA-4P

Product #:

145700

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Hardness, Total & Calcium drop count titration/EDTA test kit, range 17.1 - 342 mg/L CaCO_3 , 100 tests.

The kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

Utilises drop count titration - a simple, economical titration method

Specifications

Case Style:	D
Measuring range:	17.1 - 342 mg/L CaCO_3
Method:	Drop Count Titration/EDTA
Number of tests:	100
Parameter:	Hardness, Total & Calcium - As CaCO_3
Platform:	Drop Count Titration
Smallest increment:	17.1 mg/L
Test Kit Model:	HA-4P
What's included?:	Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.

What's included?

Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.