



## Dissolved oxygen test kit, model OX-2P

**Product #:**

**146900**

ZAR Price:  
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count titration

Replacement reagent set PN 2439802

**Utilizes drop count titration - a simple, economical titration method**

---

### Specifications

|                     |   |
|---------------------|---|
| Case Style:         | D   |
| Measuring range:    | 0.2 - 4.0 mg/L & 1 - 20 mg/L O <sub>2</sub>   |
| Method:             | Drop count titration/<br><br>Modified Winkler |
| Number of tests:    | 100   |
| Parameter:          | Dissolved Oxygen - As DO                      |
| Platform:           | Drop Count Titration                          |
| Smallest increment: | 0.2   |
| Test Kit Model:     | OX-2P   |