



## Alkalinity test kit AL-AP, 5 - 100 / 20 - 400 mg/L $\text{CaCO}_3$ 100 pcs

**Product #:**

**2444301**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists.

Method: Drop count titration

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Kit utilizes drop count titration method and contains 100 tests.

**Simple, economical titration method**

**Rugged kit contains all required reagents, solutions, and glassware**

---

### Specifications

Case Style:	B
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$ , ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Measuring range:	5 - 100 mg/L & 20 - 400 mg/L
Method:	Drop count titration/Sulfuric acid
Number of tests:	100
Parameter:	Alkalinity - (P) & Total (MO) As $\text{CaCO}_3$
Platform:	Drop Count Titration
Smallest increment:	5
	20
Test Kit Model:	AL-AP, MG-L