



Test kit, hypochlorite, model CN-HR

Product #:

2687200

ZAR Price:

Contact Hach

Ships within 3-5 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Hypochlorite Kit, Model CN-HR - Drop Count

Utilizing the simple, economical drop count titration method, this kit is used to determine high range hypochlorite solution concentrations (bleach solutions). Range: 5-15%, contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Specifications

Case Style:	D
Measuring range:	50 - 150 g/L Cl_2 (5-15%)
Method:	Drop count titration/ Thiosulfate
Number of tests:	100
Parameter:	Hypochlorite, high range - (Bleach) As Cl_2
Platform:	Drop Count Titration
Smallest increment:	0.5%
Test Kit Model:	CN-HR