



Test kit, carbon dioxide, model CA-DT

Product #:

2064100

ZAR Price:

Contact Hach

Ships within 2 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Utilizing the Hach Digital Titrator, the Model CA-DT kit is used to determine carbon dioxide concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10.1-1,000 mg/L. Contains approximately 50 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifications

| | |
|---------------------|---|
| Case Style: | F |
| Footnote: | *mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency. |
| Method: | Digital Titrator/ Sodium hydroxide Phenolphthalein |
| Number of tests: | 50 |
| Parameter: | Carbon Dioxide - As CO_2 |
| Platform: | Digital Titrator |
| Range: | 10 - 1,000 mg/L CO_2 |
| Smallest increment: | 0.1, 2.0 |
| Test Kit Model: | CA-DT |